



TideCell


Product Selection Guide

Step 1. Understand TideCell's Modules & Consumables




Step 2. Choose Your Culture Mode



Equipment



Consumables




Batch Mode




Fed-Batch Mode

Equipment




Consumables



Perfusion Mode

Step 3. Choose Your Scale based on your development status



Step 4. Choose Your Scale based on equivalency to other systems

Model	Fixed-Bed Volume (L)	Culture Surface Area (m ²)	Equivalent RB-850 (850 cm ²)	Equivalent CS-40 (25200 cm ²)	Equivalent Cytodex II (0.44 m ² /g)
BelloCell	0.1	1.65	20	0.65	3.75 g
TideCell-002	1~4	66	800	26	150 g
TideCell-010	5~10	165	2,000	65	375 g
TideCell-020	10~20	330	4,000	130	750 g
TideCell-100	50~100	1,650	20,000	650	3,750 g

Step 5. Contact us for assistance

info@cescobio.com.tw

Or local distributors